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Cable Winch Retractor acc. to Ulrich

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Infinitely variable retraction with the cable winch retractor


For precise, individual applications in the peritoneal cavity

Unique throughout the world, the cable winch technology used in the system allows surgeons to adjust retraction as desired with no preset positions when working in the peritoneal cavity. Thanks to flexible positioning for the insulated cables on the overhead bar, surgeons using the cable winch retractor can adapt the system to the situation at hand, adjusting the direction of traction and the position of the blades with exceptional precision. The retractor system is easy to mount on the upper section of the operating table, providing an optimum view of the operating site – even in difficult cases – and giving the surgeon the best possible visibility and range of motion.

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- Ideal patient care and reduced tissue trauma thanks to a wide selection of blades in a variety of shapes and sizes
- The retractor is easy to set up and position thanks to the scalable support rods and the three mounting positions for the overhead bar – suitable for operating tables of any width
- The system can be processed in a standard sterilization container thanks to the foldable design of the frame
- Considerable stability and longevity thanks to carefully selected, high-quality materials
- Insulated cables
- NEW: Now features an extra-wide overhead bar for extreme lateral cable positioning, as well as larger blades for optimum care of obese patients.

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Innovative technical improvements for model from 2018

NEW

The cable is now completely coated

- Durable and reduces costs for cable exchange
- Higher tensile load of minimum 70 kg also for obese patients

The new material of the adjusting knob is PEEK

- Avoids breakage also in case of shock of frequent sterilization

The hinge at the frame is fixed by a PEEK slider. The pin is omitted.

- No lost pins
- Easy and fast mounting of the frame

All components of the Cable Winch Retractor are easy and safe in cleaning and processing. All improved components are compatible to existing systems.

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Benefits at a glance:

Optimized access to the operation field

- The lift of the thorax provides excellent view into the operation field
- The set up of the retractor is self-explanatory and provides an access for precise surgery
- Less assistance during surgery is needed

Simple and easy in handling


- Easy processing
- Simple set-up
- Infinitely variable retraction
- Traction can be applied from an extreme lateral position

Compact system

- Only essential components to perform a liver transplantation
- High durability of the system
- Further blades and frames available

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1. The fixation element (UL7807-11) is secured to the standard rail of the OR table.

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2. The mounting support (UL7805-61/UL7805-62) is secured in the fixation element (UL7807-11).

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3. The offset blocks of the mounting support (UL7805-62) are mounted.

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4. The frame (UL7804-63 or UL7804-80) is placed on the mounting support.

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5. The cable winches (UL7801-09) are placed on the frame and provided with a blade.

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Both ends of the retractor system attach to the lateral bars of the operating table. The system consists of the following:

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**Cable winch Retractor
Complete set with blades**

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
**Cable winch Retractor extra wide
Complete set with blades**

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The system can be processed in a standard container and tray.

GEISTER®

Seilzugsystem
Cable Winch System

Länge length	Breite width	Tiefe depth	Artikel-Nr. Product No.
88	63	69	UL7801-09

Feinstufig justierbar
Infinitely variable adjustment

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Reparatur-Set
Repair-Kit

Artikel-Nr. Product No.
UL7801-63

Inbusschlüssel und Loctite 5g
Hex-wrench and Loctite 5g

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GEISTER Medizintechnik GmbH