







ValveGate™*PRO*

Efficient surgery. Less invasive. Optimum results.

Achieve smaller approaches. Handle greater challenges.

Minimally invasive approaches, performed through small incisions, are becoming state-of-the-art in an effort by healthcare professionals to provide each patient with the best possible treatment. The most common types of minimally invasive cardiac surgery (MICS) are valve surgeries, including repairs and replacements, and even more challenging, coronary artery bypass. Such smaller incisions without a full sternotomy are expected to leave smaller scars with less trauma and bleeding.

Celebrating 20 years of excellence in MICS.

20 years experience and close cooperation with opinion leaders in Minimally Invasive Cardiac Surgery resulted in the development of our latest and most advanced generation of ValveGate™ MICS instruments... the all new ValveGate™PRO line!

Features:

- Ultra-lightweight and perfect balance
- Completely made from surgical stainless steel without any plastic or aluminium components
- Full injection straight-flush cleaning concept
- Adjustable jaw and handle opening according users preference for less fatigue and time savings
- Scissors designed to cut even through hardest calcified tissues
- Needle Holders firmly grasp and drive even the largest needles
- Precise anastomosis suturing by finger turning and not by wrist bending
- New high precision ratchet
- Soft touch handles with elastic springs
- Ergonomic FAIRGRIP™ handles with large dimples
- Solid construction and optimized surfaces for maximum durability, wear and corrosion resistance
- Color coded handles (Gold= Needle Holders, Black= Scissors, Silver= Graspers, all others).





ValveGate[™] PRO





ValveGate[™] *PRO* Scissors

15° curved

V 38-283U

А	30-7030	15 Curveu
В	38-7831	Jameson Scissors, curved,
		tenotomy tips
С	38-7832	20° curved, single action,
D	38-7833	robust, 30° curved, 7 mm shaft
Ε	38-7835	30° curved
F	38-7836	40° curved, single action
G	38-7840	70° curved
Н	38-7842	Ligature Scissors, 70° curved wave-cu
1	38-7860	Potts Scissors 125°, 3.4 mm shaft
J	38-7862	Potts Scissors 45°, 3.4 mm shaft
K	38-7863	Potts Scissors 125°, 5 mm shaft
L	38-7861	Potts Scissors 30°, heavy pattern

ValveGate[™] *PRO* Clip Applying Forceps

M	38-7460	small, for use with SLS/Vitalitec* Clips
	38-7462	medium, for use with SLS/Vitalitec* Clips

5 mm shaft diameter, unless otherwise specified. Standard working length 25 cm/10" for minimally invasive valve surgery. For different length add "B", "C", "S", "L" or "XL" behind article number.

B = 12 cm/5" C = 15 cm/6" S = 17 cm/6 3 /4" L = 30 cm/12" XL = 35 cm/14"

^{*}Trademarks property of their respective trademark owner.



ValveGate[™] PRO





ValveGate™ PRO



ValveGate™ *PRO* Needle Holders with ratchet

А	30-1199	10 curvea, Grile-wood
В	38-7800	TC straight
С	38-7802	TC curved
D	38-7803	TC curved, short Colvin jaw
Ε	38-7804	TC curved, del., Micro Crile-Wood
		jaw for 5-0/6-0 needles
F	38-7806	Mini TC curved, 2.7 mm shaft
G	38-7805	TC Ryder straight
Н	38-7807	TC angled up

ValveGate[™] *PRO* Dissecting Forceps

1	38-7823	Leaflet Grasper, long 28 mm,
		with ratchet
1	38-7823.NC	Leaflet Grasper, long 28 mm,
		without ratchet
J	38-7828	Gemini Dissecting Forceps,
		smooth dilating tips 25 mm

ValveGate™ PRO Knot Pusher

K	38-7850	push to close
L	38-7851	push to open

ValveGate™ PRO Mini IMA Holding Forceps

M 38-7878 Jaw 0,8x7 mm, upwards, 1x2 teeth, Sapphire™, 3,4 mm shaft

5 mm shaft diameter, unless otherwise specified. Standard working length 25 cm/10" for minimally invasive valve surgery. For different length add "B", "C", "S", "L" or "XL" behind article number.

B = 12 cm/5" C = 15 cm/6" S = 17 cm/6 3 4" L = 30 cm/12" XL = 35 cm/14"



Wet-Set for transport and storage for ValveGate™ Instruments

Keeps the distal end moist and facilitates the reprocessing. Protects the instrument tip during transport.

00-2250 Set: 3 tubes and caps (with hole) 00-2251 Set: 20 tubes and caps (with hole)

00-2252 Set: 20 tubes and caps (without hole)


00-2253 Set: 20 tubes and 10 caps without

hole and 10 caps with hole



ValveGate[™] PRO





Α	38-7935	10° angled
В	38-7936	25° angled
	38-7934	straight

ValveGate™ *PRO* Graspers C 38-7824 Resano, strong, jaw for valves, straight

D	38-7929	Russian, jaw 6,8x5,3 mm
F	38-7567	Lymph Node Grasper angled up, jaw 7 mm ovale, atraumat

ValveGate™ PRO Gentle Graspers

F	38-7816	extra wide jaw, straight, TC
G	38-7826	extra wide jaw, straight, 2x3 teeth,
		4,8x14 mm, with extra fine teeth

ValveGate™ PRO Graspers with DeBakey Jaws, 1x2 teeth

Н	38-7820	narrow, straight, 1.5x11 mm
1	38-7821	narrow, angled upwards, 1.5x14 mm
J	38-7822	narrow, angled sidewards, 1.5x18 mm, vertically opening
K	38-7825	wide jaw, 2.8x11 mm, straight
L	38-7827	delicate, 1.5x11 mm, 3 mm malleable shaft
M	38-7881	narrow, 1.5x15 mm, straight
Ν	38-7917	narrow, 1.5x15 mm, slightly curved
0	38-7918	narrow, 1.5x15 mm, angled sidewards 10°
Р	38-7919	narrow, 1.5x15 mm, angled sidewards 25°
Q	38-7922	narrow, angled sidewards, 1.5x18 mm, horizontally opening
R	38-7928	long, 2,8x25 mm teeth, 40 mm jaw

5 mm shaft diameter, unless otherwise specified. Standard working length 25 cm/10" for minimally invasive valve surgery. For different length add "B", "C", "S", "L" or "XL" behind article number.

B = 12 cm/5" C = 15 cm/6" S = 17 cm/6%" L = 30 cm/12" XL = 35 cm/14"





ValveGate™





ValveGate[™] Mohr Suture Ruler

(Chordae Tendinae Measuring and Knotting

Device for ePTFE sutures)

03-5409 ValveGate™ Suture Ruler

3-60 mm, 35 cm/13¾"

ValveGate[™] Mohr Knotting Device

for chordal reconstruction with ePTFE sutures 03-5420

Shibata Chordae System

knotting device for facilitating chordal loop making

03-5425 15-26 mm 03-5425L 11-30 mm

Adams-Yozu Suture Gauge

03-5432

Adams-Yozu Knotting Caliper

with weight holder 03-5430



ValveGate[™] Accessories



34-7471







Α	34-7486	ValveGate™ Fine Nerve and Vessel Hook, ball point, 90°	35 cm/13¾"
	34-7487	ValveGate™ Fine Nerve and Vessel Hook, ball point, 70°	35 cm/13¾"
В	34-7485	ValveGate™ Nerve and Vessel Hook	35 cm/13¾"
С	34-7484	ValveGate™ Vessel Hook	35 cm/13¾"
D	34-7475	ValveGate™ Suture Catcher fine, 1.8 mm	35.5 cm/14"
	34-7476	ValveGate™ Suture Catcher extra fine, 1.5mm	35.5 cm/14"
Ε	34-7480	ValveGate™ Magnetic Retriever	35 cm/13¾"
F	34-7471	ValveGate™ Holder for Scalpel blades #11	35 cm/13¾"
		(without blades)	
	16-2011	Scalpel Blades #11 (pack of 100, sterile)	
G	34-7470	ValveGate™ Holder for Mini Blades (without blades)	37 cm/141/2"
Н	34-7474	ValveGate™ Protected Mini Blade Holder	37 cm/141/2"
	16-2065	Mini Scalpel Blades #65 (pack of 10, sterile)	
- 1	34-7465	ValveGate™ Valve Pusher	35 cm/13¾"
J	34-7494	ValveGate™ Knot Pusher, both sided	35 cm/13¾"
K	34-7493	ValveGate™ Knot Pusher with blue silicone handle	31.5 cm/121/2"
K	34-7495	ValveGate™ Knot Pusher	35 cm/13¾"
L	34-7496	ValveGate™ Knot Pusher with disc	35 cm/13¾"



GEISTER®

ValveGate™ Endoscopes, Trocars, Sleeves

Endoscopes				
Art.No.	Diameter	Direction		
	Working length	of view:		
34-0960.17	Ø 5.0 mm, 17,5 cm	0°		
34-0961.17	Ø 5.0 mm, 17,5 cm	30°		
34-0961.23	Ø 5.0 mm, 23 cm	30°		
34-0960.AC	Ø 5.0 mm, 30 cm	0°		
34-0961.AC	Ø 5.0 mm, 30 cm	30°		
34-0962.18	Ø 10.0 mm, 18 cm	0°		
34-0967.18	Ø 10.0 mm, 18 cm	30°		
34-0962.23	Ø 10.0 mm, 23 cm	0°		
34-0967.23	Ø 10.0 mm, 18 cm	30°		
34-0962.AC	Ø 10.0 mm, 31,5 cm	0°		
34-0967.AC	Ø 10.0 mm, 31,5 cm	30°		
34-0962.35	Ø 10.0 mm, 35 cm	0°		
34-0967.35	Ø 10.0 mm, 35 cm	30°		
34-0962.43	Ø 10.0 mm, 43 cm	0°		
34-0967.43	Ø 10.0 mm, 43 cm	30°		
ACMI-, Wolf- and Storz-lightport adapter, standard				

eye-piece, autoclavable, S/N, single packed.

ValveGate™ Flushing Needles

Art.No.	Diameter	Working length
18-0050	Ø 2.0 mm	305 mm / 12"
18-0050L	Ø 2.0 mm	340 mm / 13½"
18-0051	Ø 3.0 mm	305 mm / 12"

Suction Tips Art.No.	Description	Length of Basket	Tube i.d.	Size
18-4530	Vent - Bullet Heavy	40 mm	³ / ₁₆ "	
18-4540	Vent - Bullet Double Lumen	42 mm	1/4"	
18-4520	Vent - Ross	33 mm	³ / ₁₆ "	Adult
18-4521	Vent - Ross	23 mm	³ /18"	Child
18-4522	Vent - Ross	15 mm	1/8"	Infant

Thoracic Trocars with Sleeves					
Art.No.		Diameter, Working length			
34-0855	rigid	Ø 6 mm, 50 mm			
34-0852	rigid	Ø 11 mm, 50 mm			
34-0851	rigid	Ø 13 mm, 50 mm			
34-7510	flexible	Ø 6 mm, 70 mm			
34-7510.ST	flexible Spare Tubes, Pack of 5 pcs.	Ø 6 mm, 70 mm			
34-7511	flexible	Ø 6 mm, 100 mm			
34-7511.ST	flexible Spare Tubes, Pack of 5 pcs.	Ø 6 mm, 100 mm			
34-7500	flexible	Ø 11 mm, 100 mm			
34-7500.ST	flexible Spare Tubes, Pack of 5 pcs.	Ø 11 mm, 100 mm			
34-7501	flexible	Ø 11 mm, 70 mm			
34-7501.ST	flexible Spare Tubes, Pack of 5 pcs.	Ø 11 mm, 70 mm			
34-0860	Adams-Yozu Introducer, blunt	Ø 2.3 mm			
Trocars for	Indroducer Sleeves, metal, smooth				
pyramidal po	pint	Diameter			
34-0956		Ø 5,5 mm			
34-0903		Ø 10 mm			
34-0923		Ø 11 mm			
34-0925		Ø 13 mm			
conical point	1				
34-0957		Ø 5,5 mm			
34-0924		Ø 11 mm			


Indroducer Sleeves

Art.No.		Diameter, Working length
34-0952	Metal shaft, smooth	Ø 5,5 mm, 103 mm
34-0952.PL	Plastic shaft, smooth	Ø 5,5 mm, 92 mm
34-0952.WC	White seal for Ø 5,5 mm	
34-0952.VV	Silicone valve, pack of 5 pieces, for Ø 5,5 mm	
34-0951	Metal shaft smooth	Ø 11 mm, 105 mm
34-0951.PL	Plastic shaft, smooth	Ø 11 mm, 97 mm
34-0951.WC	Blue seal for Ø 11 mm	
34-0950	Metal shaft smooth	Ø 13 mm, 102 mm
34-0950.PL	Plastic shaft, smooth	Ø 13 mm, 94 mm
34-0950.WC	Green seal for Ø 13 mm	
34-0951.VV	Silicone valve, pack of 5 pieces, for Ø 11 and Ø 13	3 mm



ValveGate[™] Accessories







ValveGate™ Transthoracic Occlusion Clamp with DeBakey Jaws, with ratchet, sliding shaft

34-7301	100 mm (1x2) regular jaw, working length 29 cm	36 cm
34-7302	50 mm (1x2) regular jaw, working length 22 cm	32 cm
34-7303	75 mm (1x2) regular jaw, working length 24,5 cm	34 cm
34-7306	100 mm (2x3) wide jaw, working length 26,5 cm	36 cm
34-7311	100 mm (2x3) wide jaw, acutely curved, working length 26,5 cm	35 cm

	_	•		
Α		40-2370.04	4mm, straight	20cm
В		40-2371.04	4mm, 50° cutting upward	20cm
С		40-2371.04D	4mm, 50° cutting down	20cm
Α		40-2372.04	4mm, straight	30cm
В		40-2373.04	4mm, 50° cutting upward	30cm
С		40-2373.04D	4mm, 50° cutting down	30cm


ValveGate™ Valve Rongeur/Biopsy Forceps

34-7450	3mm, straight	30 cm
34-7450.Z	3mm, straight, detachable	30 cm
34-7451	3mm, 30° upwards	30 cm
34-7451.Z	3mm, 30° upwards, detachable	30 cm

ValveGate™ Bulldog Clamps and Applying Forceps

20-0351	Bulldog Clamp with DeBakey Jaws, straight	9 cm/31/2"
20-0352	Bulldog Clamp with DeBakey Jaws, straight	10.5 cm/4"
20-0378	Bulldog Clamp Applying Forceps, jaw 45° angled up	24 cm/91/2"
34-7280	Endoscopic Bulldog Clamp Applying Forceps, angled	31 cm/121/4"
34-7281	Endoscopic Bulldog Clamp Applying Forceps, straight	31 cm/121/4"

Release add "K" behind catalog reference if black ceramic coating is desired.



ValveGate™ FlexCrown™








FlexCrown™

29-1372 FlexCrown[™] Retractor Set in tray consisting of:
29-1375 FlexCrown[™] Mitral Spring Retractor*
29-1375.S FlexCrown[™] Mitral Spring Retractor, small*
29-1373 FlexCrown[™] Suture Catcher
29-1374 FlexCrown[™] Introduction Forceps
77-3645.FC Medium multisystem-basket container

with mounts for FlexCrown™ with lid,

240 x 240 x 45 mm







ValveGate[™] Holders





ValveGate™ Set of 6 Tricuspidal Blades

29-1379	consi	sting of	
ArtNo.	qty.		length x width mn
29-1379.1	1	Tricuspidal Blade with slot	25,5 x 19
29-1379.2	1	Tricuspidal Blade with slot	25,5 x 24,5
29-1379.3	1	Tricuspidal Blade with slot	25,5 x 38
29-1379.4	1	Tricuspidal Blade with slot	38,5 x 34,5
29-1379.5	1	Tricuspidal Blade with slot	51 x 25,5
29-1379.6	1	Tricuspidal Blade with slot	25,5 x 48
29-1379.STK	1	Holding Road short	
29-1379.STL	1	Holding Road long	
29-1379.HZ	1	Introducing Forceps	

ValveGate™ Set of 6 Mitral Blades

29-1380	consi	sting of	
ArtNo.	qty.		length x width mm
29-1380.1	1	Mitral Blade with slot	40 x 35
29-1380.2	1	Mitral Blade with slot	60 x 35
29-1380.3	1	Mitral Blade with slot	50 x 35
29-1380.4	1	Mitral Blade with slot	50 x 45
29-1380.5	1	Mitral Blade with slot	60 x 45
29-1380.6	1	Mitral Blade with slot	70 x 45
29-1380.STK	1	Holding Road short	
29-1380.STL	1	Holding Road long	

Introducing Forceps

Hand-Held Deep Atrial Retractor

29-1380.HZ

30-0409 working length 240 mm, 25 mm width



Iron-Assistant™ Pneumatic Endoscope Holder and Trays





29-5000 Iron-Assistant™ Pneumatic Endoscope Holder

for direct connection with the compressed air supply in the OR (working pressure 5 bar up to 7 bar). Enables instantaneous and unlimited one-hand manipulation of an endoscope or tube shaft instrument within the operating field followed by total steadiness and stability under load once in the "locked" position. Easy handling and highest reliability due to the one-hand concept with only 3 hinges and the use of titanium components for weight reduction. The pneumatic arm can be easily attached to the OR table rail and is autoclavable completely or partially (if draping is preferred) for quick sterile set-up.





Large special Trays for Iron-Assistant™

77-3623 536 x 240 x 85 mm, with lid,

for 29-1481, 29-1481/I, 29-1482, 29-1482/I

77-6162.01 535 x 340 x 95 mm for 29-1481G

77-6162.01LD Lid for 77-6162.01

77-6164 644 x 249 x 140 mm, with lid, for Pneumatic 29-5000



Iron Assistant™ Endoscope and Instrument Holders





Iron Assistant™

29-1481G	Instrument Holder, extra long straight bar
29-1481	Instrument Holder, L-shaped bar
29-1481/I	ditto, insulated for use with HF-surgery
29-1482	Endoscope Holder, straight bar
29-1482/I	ditto, insulated for use with HF-surgery
29-1482L	Endoscope Holder L-shaped bar
29-1482L/I	ditto, insulated for use with HF-surgery
29-1483	Instrument Holder, L-shaped bar, double arm
	(Left lateral LVAD implantation)

Iron Assistant™ with quick connect

29-1800	Holder, 360° rotatable
29-1804	Table clamp, heavy duty
29-1808	Instrument Clamp
29-1809	Endscope Clamp

Suitable Trays see page 25



$ValveGate^{\tiny{\mathsf{TM}}}\ Mini\text{-}Thoracotomy\ Retractor}$





29	\circ	\sim	\sim
Ju		ч	n



Retractor complete including turnable crossbar (detachable) for accessories						
ArtNo.	description		Α	В	С	D
29-1395	straight arms, tubular blades		60 mm	40 mm	90 mm	120 mm
29-1390	curved arms, tubular blades, with stro	ong spreading bar	60 mm	40 mm	120 mm	75 mm
29-1390.45	curved arms, tubular blades, with stro	ong spreading bar	45 mm	40 mm	90 mm	75 mm
29-1390.80	curved arms, tubular blades, with stro	ong spreading bar	80 mm	40 mm	120 mm	75 mm
with hinged a	arms, exchangeable blades	Blade Pair				
29-1396	with 2 x 2 parallel blades*	29-1396.40	35 mm	35 mm	165 mm	100 mm
		29-1396.50	50 mm	35 mm		
29-1396E2	with 2 x 2 tubular blades**	29-1396E.50	50 mm	40 mm	165 mm	100 mm
		29-1396E.70	70 mm	40 mm		
29-1396E3	with 3 x 2 tubular blades**	29-1396E.30	30 mm	40 mm	165 mm	100 mm
		29-1396E.50	50 mm	40 mm		
		29-1396E.70	70 mm	40 mm		
A	December ded Value Cata IM					

Accessories - neconinienceu valvedate		
ArtNo.	description	length, width
29-1397.06	Suture Ring with removable springs, detachable	Ø 6 mm
29-1397.08	Suture Ring with removable springs, detachable	Ø 8 mm
29-1397.06FH	Spring Stretchers (2/set) for Suture Rings 29-1397.06 and	129-1397.08
29-1395.TB	Extra turnable crossbar	
Accessories -	Atrial Hooks with 6 mm Ø rod	lenath, width

Accessories -	Atrial Hooks with 6 mm Ø rod	length, width
29-1500.BH	malleable	25 cm, 20 mm
29-1500.HS	Cooley, rigid, curved basket	23 cm, 20 mm
29-1501.HS	Cooley, rigid, straight basket	25 cm, 20 mm
29-1500.HS/D	Cooley, rigid, curved basket, long shaft	31cm, 20 mm
29-1428	Fixation Clamp for atrial hooks	Ø 6 mm
29-1437	Multi-Purpose Spatula	29 cm, 20 mm
10-8000	Blade Changing Forceps (is including in 29-1396, 29-1396)	96E2, 29-1396E3)

Additional Acc	essories	iength, wiath
30-0740	Cone Retractor, deep pattern	14,5 cm/5 3/4"
30-0546.16	Memory malleable Spatula	18 cm/7", 16 mm
30-0546.19	Memory malleable Spatula	18 cm/7", 19 mm
30-0546.24	Memory malleable Spatula	18 cm/7", 24 mm
30-0547.20	Memory malleable Spatula	20 cm/8", 20 mm

^{*} shown with 29-1428 and 29-1500.HS (not included) ** also available with fenestrated blades. Add "M" behind article number.

ValveGate™ Retractors













29-0938



Adams-Yozu Retractor					
ArtNo.	description	Α	В	С	D
29-0930	fix blade	55 mm	50 mm	100 mm	95 mm
	exchangeable blade	65 mm	50 mm		
smal version					
29-0929	fix blade		40 mm	100 mm	95 mm

Accessories Art.-No.

description

29-0933	Adams-Yozu Mini Arm with M6 screw connector
29-0934	Adams-Yozu Mini Arm with quick connector for spreading bar
29-0935	Attachment and Extension blade

Depressor b	total le	total length			
29-0937.18	XXS, Pediatrie	46 cm	85 mm	18 mm	
29-0937.23	XS	44 cm	70 mm	23 mm	
29-0937.28	S	44 cm	70 mm	28 mm	
29-0937.35	M	44 cm	70 mm	35 mm	
29-0937.45	L	44 cm	70 mm	45 mm	
29-037.45L	XL	47 cm	90 mm	45 mm	
Spatula malleable					
29-0938		36 cm	70 mm	20 mm	



ValveGate[™] Soft Tissue Protectors




ValveGate[™] Soft Tissue Protector

The ValveGate™ Soft Tissue Protector for Minimally Invasive Valve Surgery is a sterile, single patient use wall access device consisting of 2 overlaid plastic rings that are interconnected by means of a silicone rubber membrane. The bottom ring is a flexible ring made with a shape memory alloy to be inserted in situ.

The special features are:

- clear and unimpeded access to surgical site in minimally invasive valve or bypass surgery, MIDCAB, V.A.T.S., transcatheter or transapical valve replacement a.o.
- protects wound edge and surrounding tissue
- usable with or without retractor
- can be used in combination with incision foil
- reduces risk of embolic fat particles on instruments, sutures and knot pushers
- improves visibility and exposure
- self-holding, for incisions from 6 cm to 9 cm (regular size) and from 3 cm to 5 cm (small size)
- very easy to insert in situ
- high flexibility
- latex-free silicone rubber

Art.-No. 29-1399 regular size 120 x 100 mm
Art.-No. 29-1399M medium size 100 x 100 mm
Art.-No. 29-1399S small size 70 x 70 mm









For easy placement simply push the flexible ring downwards as illustrated and insert it into the wound cavity. Then push the reversed portion carefully down as well and expand it. Remove it in the opposite way by inserting some fingers, grabbing a part of the flexible ring, and pulling it out carefully. Dispose after use.







ValveGate™ Valve Sizer






ValveGate™ Valve Sizer

03-5465.00 ArtNo.	Valve Sizers complete consisting of : description
03-5465	Valve Sizer holder without sizers
03-5465.25	Valve Sizer #25
03-5465.27	Valve sizer #27
03-5465.29	Valve Sizer #29
03-5465.31	Valve Sizer #31
03-5465.33	Valve Sizer #33
77-3641	Mini multisystem-basket container with lid,
	with variable mount, 126 x 14 x 39 mm









System Trays











77-3620 ValveGate[™] PRO System Tray basket container with filter

Special Tray DIN with filter to connect 8 instruments

480 x 240 x 45 mm, with lid

77-3621 ValveGate[™] *PRO* System Tray basket container

Special Tray DIN with variable bracket and mounts

480 x 240 x 45 mm, with lid

77-3600 Plastic Speciality Containers

Convertible ValveGate[™] & Laparoscopy, Base w. mat, removable insert, 478 x 203 x 130 mm, clear lid

77-3623 ValveGate[™] *PRO* System Tray

large special Tray for Iron-Assistant™ (29-1481, 29-1481/I, 29-1482, 29-1482/I)

536 x 240 x 85 mm, with lid

00-3800 Charly ValveGate Trainer

with USB-camera and -light

Suitable Boxes and Accessories see special catalog "Trays and Containers"





ValveGate™ classic line **Scissors** Fig. Art.-No. Fig. Art.-No. Fig. Art.-No. Fig. Art.-No. Fig. Art.-No. Fig. Art.-No. 34-7830 34-7832 34-7835 34-7840 34-7860 11 34-7958 34-7833 34-7842 10 34-7862 34-7831 34-7836 12 34-7861 Fig. Art.-No. 13 34-7816 21 34-7875S 34-7919 15 34-7821 17 34-7824 19 34-7826 23 34-7877S 25 34-7917 27 14 34-7820 16 34-7822 18 34-7825 20 34-7827 22 34-7876S 24 34-7881 26 34-7918 28 34-7929 **Needle Holders** Fig. Art.-No. Fig. Art.-No. Fig. Art.-No. Fig. Art.-No. Fig. Art.-No. 29 34-7799 34-7802 33 34-7804 34-7805 34-7808S 30 34-7800 34-7803 34 34-7806 36 34-7807 34-7952 **Dissecting Forceps** Fig. Art.-No. Fig. Art.-No. Fig. Art.-No. Fig. Art.-No. 39 34-7823 41 34-7878 42 34-7460 44 34-7850 40 34-7828 43 34-7462 45 34-7851



REMARKABLE

Black and White Special Edition

38-7800/WK A variety of models is also optionally available in a new design of Black & White Special Edition (extra charge applies), please add "/WK" behind the item-number for ordering.

Pandor[™] Special Edition

38-7820.DF

A variety of models is also optionally available in the attractive PANDOR™ Black & Gold Edition (extra charge applies), please add ".DF" behind the item-number for ordering.

GEISTER Medizintechnik GmbH | Föhrenstrasse 2 | 78532 Tuttlingen/Germany Tel +49 74 61/966 240 | Fax +49 74 61/966 24 22 | Mail info@geister.com | www.geister.com



